LabX Balance Software

Digitalize Your Weighing Processes



Seamless Integrated Dataflow

LabX integrates fully with laboratory information systems, such as LIMS and chromatography software, and sends your weighing results directly to where the weight values are required for subsequent analyses.



Automatic Documentation

LabX automatically saves weighing data, including all metadata, in a secure, centralized database at the time of creation. Customized result reports can be generated at any time, including automatically at the end of a procedure.




Central Balance Management

User profiles, weighing tasks and balance configurations can be managed centrally. LabX monitors each balance's test and calibration history and automatically initiates routine tests or blocks the balance when necessary.



Improved Efficiency

Thanks to the intuitive instructions on the balance display and automatic documentation, using LabX saves a considerable amount of time and significantly reduces errors. Efficiency is improved and output increases.



Guided Digital Workflows

And Automatic Documentation

Start your weighing task with just one touch on your balance display. LabX™ Balance guides you step-by-step through your SOP with clear onscreen instructions at the balance. All results and process metadata are documented and saved in a centralized database, calculations are performed automatically and transcription errors are eliminated. LabX provides a seamless flow of data throughout the entire analysis workflow, guaranteeing a high degree of process security. LabX can also support you with the fulfilment of the FDA ALCOA+ requirements for data integrity.

Welcome to the digital laboratory where your data is available from anywhere 24/7 and you are always audit-ready.



One System

for Multiple Instruments

LabX software enables up to 60 METTLER TOLEDO instruments to be fully connected via your internal network. Enjoy seamless operations and efficient workflows with Automatic Balances, Balances, Titrators, pH Meters, UV Vis Spectrophotometers, Density Meters, Refractometers and Melting Point Instruments.



LabX Balance Starter Packs

LabX Balance Express

- Installation on a single PC
- Connect up to 3 balances

LabX Balance Server

- Installation on multiple PCs
- Can be connected in a distributed network
- Connect up to 60 balances

Starter Pack LabX Balance Express 1 Automatic Balance / Balance Instrument License

Order no. 11153120

1 Automatic Balance / Balance	
Instrument License	
User Management	
Auto Import / Export	

Add Optional Licenses User Management 11153102 Statistica Auto Import / Export 11153105 Report De Regulation* 11153104 System In Product Database 11153103 STAREX In

Statistical Evaluation	30262374
Report Designer	11153106
System Integration	11153107
STAReX Integration	30486416

Add Additional Instrument Licenses	
Balance / Balance Instrument License	11153220

^{*} Also includes the option User Management

LabX Software Services

METTLER TOLEDO offers a full range of services for LabX Balance Software:

- Installation
- Training
- Validation
- Integration
- Workflow Development

Complete your system with a Software Care contract to stay up to date with all the latest innovations, benefit from consistent performance and ensure regulatory compliance. Software Care provides you immediate access to new releases and professional support.

For further details, please refer to the LabX Services brochure or contact your local METTLER TOLEDO representative.

METTLER TOLEDO Group

Laboratory Weighing Local contact: www.mt.com/contacts

Subject to technical changes
© 12/2020 METTLER TOLEDO. All rights reserved
30207869B
Group MarCom RITM692067 LS

www.mt.com

For more information