

Mettler-Toledo GmbH
Laboratory & Weighing Technologies

Address Im Langacher, CH-8606 Greifensee, Switzerland
Phone +41-44-944 36 36
Fax +41-44-944 31 90
E-mail support.labtec@mt.com

Release date June 22nd, 2017

Software Release Chart - SerialPortToKeyboard

Date	Software Version	Issue Description
22.06.2017	1.11	Compatibility with MS Office 365™ and Office 2016™
01.05.2015	1.10	Compatibility with MS Excel on Windows 10™
18.04.2013	1.00	Initial version

In each new software version are bug fixes included. The Information contained in this document has been collected to the best of our knowledge and represents the latest update. In addition, Mettler-Toledo GmbH will in no event be liable for consequences of any kind arising out of, or in connection with, the use of this SW or this document.

Operating System Compatibility for Software Version 1.11

Software Description	Version	Service Pack
Operating System (32/64 bit)	Microsoft® Windows 7™	SP1
Operating System (32/64 bit)	Microsoft® Windows 8™	SP1
Operating System (32/64 bit)	Microsoft® Windows 10™	
Operating System (32/64 bit)	Microsoft® Windows 11™	

Notes:

Microsoft® Windows Operating Systems 95 / 98 / ME / NT / 2000 / XP / Vista are no longer supported.

Installation Note

Software must be installed with a supervisory-level login user name, such as "Administrator". Right-click the installer .EXE and select 'Run as Administrator'.

Executing the installation file a second time and confirming the 'Confirm Uninstall' dialog leads to the complete removal of all files installed except the configuration files.

Company name: Mettler-Toledo GmbH
Laboratory & Weighing Technologies

Date: June 22nd, 2017

Function: Product Manager SPG Balances

Signature: Michael Jucker

